

# TOPIVAC NEGATIVE PRESSURE TOPICAL VACUUM WOUND CARE THERAPY DEVICES

## Technical Specifications / Datasheet



MODELS	HAND T-NPWT CLASSIC	HAND T-NPWT DIGITAL IRRIGATION	MEDIUM V2	MEDIUM CLINIC / V4
TYPES	Digital Hand Vacuum Assisted(NPWT) Wound Closure Therapy Device	Digital Hand Vacuum Assisted Irrigation Wound Closure Therapy Device	Medium V2 Professional Type Hybrid Topical Wound Treatment Device	Medium Clinic / V4 Professional Type Hybrid Topical Wound Treatment Device
CLASS	MDD: Class IIa	MDD: Class IIa	MDD: Class IIa	MDD: Class IIa
SIZE (mm)	Width:150 Length:190 Height:200	Width:150 Length:190 Height:200	Width:150 Lengt:190 Height:200	Width:150 Lengt:190 Height:200
SCREEN	LCD 5" Touch Screen	LCD 5" Touch Screen	LCD 5" Touch Screen	LCD 5" Touch Screen
POWER	2,4W (stand-by in use)	2,4W (stand-by in use)	3 W (stand-by in use)	3 W (stand-by in use)
VOLTAGE	220 V - 50/60 Hz / 13.8 VDC	220 V - 50/60 Hz / 13.8 VDC	220 V - 50/60 Hz / 13.8 VDC	220 V - 50/60 Hz / 13.8 VDC
ADAPTER	Input 220 VAC 50-60 Hz output 13.8 VDC	Input 220 VAC 50-60 Hz output 13.8 VDC	Input 220 VAC 50-60 Hz output 13.8 VDC	Input 220 VAC 50-60 Hz output 13.8 VDC
WEIGHT	1,650 Kg.	1,950 Kg.	2,050 Kg.	2,100 Kg.
VACUUM	-20 / -200 mmHg	-20 / -200 mmHg	-20 / -200 mmHg	-20 / -200 mmHg
VOLUME	55 dB, ±3	55 dB, ±3	55 dB, ±3	55 dB, ±3
ALARMLAR	Visual & Audible Warning	Visual & Audible Warning	Visual & Audible Warning	Visual & Audible Warning
FUNCTIONAL SPECIFICATION	<ul style="list-style-type: none"> <li>• Intermittent and Incremental Vacuum Therapy</li> <li>• Classic Negative Pressure Vacuum Therapy</li> </ul>	<ul style="list-style-type: none"> <li>• Intermittent and Incremental Vacuum Therapy</li> <li>• Classic Negative Pressure Vacuum Therapy</li> <li>• Moistening</li> <li>• Irrigation with Bidistilled Water</li> <li>• Special Irrigation</li> </ul>	<ul style="list-style-type: none"> <li>• Intermittent and Incremental Vacuum Therapy</li> <li>• Classic Negative Pressure Vacuum Therapy</li> <li>• Moistening</li> <li>• Irrigation with Bidistilled Water</li> <li>• Special Irrigation</li> </ul>	<ul style="list-style-type: none"> <li>• Intermittent and Incremental Vacuum Therapy</li> <li>• Classic Negative Pressure Vacuum Therapy</li> <li>• Moistening</li> <li>• Irrigation with Bidistilled Water</li> <li>• Special Irrigation</li> </ul>

Compatible Topical Vacuum Treatment Dressing Sets for Devices: Multidress, Multicase, Canister, Ozonbag, Topiset